

COLORODENT BASIC B3

COLORI PER ZIRCONIA AD USO DENTALE
ZIRCONIA COLORS FOR DENTISTRY
ISO 6872 STANDARD

CONTIENE: DIETILEN GLICOL
IT CONTAINS: DIETILENE GLICOLE

COMPOSIZIONE %:
COMPOSITION %

TC100901 3X	TC100903 3X	TC100904 3X	TC100905 3X	TC100999
38%	8%	8%	9%	37%
0.32% Pr + 66% H ₂ O	0.09% Fe-0.09% Cr + 66% H ₂ O	0.25% Co+ 66% H ₂ O	0.10% Cr-1.02% Al + 66% H ₂ O	C ₄ H ₁₀ O ₃

AVVERTENZE
WARNINGS
Xn



NOCIVO

NOCIVO PER INGESTIONE
R22
HARMFUL TO RECORDING