

COLORODENT BASIC D2

COLORI PER ZIRCONIA AD USO DENTALE
ZIRCONIA COLORS FOR DENTISTRY
ISO 6872 STANDARD

CONTIENE: DIETILEN GLICOL
IT CONTAINS: DIETILENE GLICOLE

COMPOSIZIONE %:
COMPOSITION %

TC100901 3X	TC100903 3X	TC100904 3X	TC100905 3X	TC100999
12,5%	14,5%	4,5%	4,5%	64%
0.32%Pr + 66% H ₂ O	0.09%Fe-0.09%Cr + 66% H ₂ O	0.25% Co+ 66% H ₂ O	0.10%Cr-1.02%Al + 66% H ₂ O	C ₄ H ₁₀ O ₃

AVVERTENZE
WARNINGS
Xn



NOCIVO

NOCIVO PER INGESTIONE
R22
HARMFUL TO RECORDING