

# COLORODENT BASIC D4

**COLORI PER ZIRCONIA AD USO DENTALE  
ZIRCONIA COLORS FOR DENTISTRY  
ISO 6872 STANDARD**

**CONTIENE: DIETILEN GLICOL  
IT CONTAINS: DIETILENE GLICOLE**

**COMPOSIZIONE %:  
COMPOSITION %**

<b>TC100901 3X</b>	<b>TC100903 3X</b>	<b>TC100904 3X</b>	<b>TC100905 3X</b>	<b>TC100999</b>
<b>48%</b>	<b>10%</b>	<b>10%</b>	<b>15%</b>	<b>17%</b>
0.32% Pr + 66% H <sub>2</sub> O	0.09% Fe-0.09% Cr + 66% H <sub>2</sub> O	0.25% Co+ 66% H <sub>2</sub> O	0.10% Cr-1.02% Al + 66% H <sub>2</sub> O	C <sub>4</sub> H <sub>10</sub> O <sub>3</sub>

**AVVERTENZE  
WARNINGS**  
Xn



**NOCIVO**

**NOCIVO PER INGESTIONE  
R22  
HARMFUL TO RECORDING**